

Contents lists available at ScienceDirect

Journal of Virological Methods



journal homepage: www.elsevier.com/locate/jviromet

Protocols

Sensitive detection of African swine fever virus using real-time PCR with a 5' conjugated minor groove binder probe

John McKillen^{a,*}, Michael McMenamy^b, Bernt Hjertner^b, Francis McNeilly^a, Åse Uttenthal^c, Carmina Gallardo^d, Brian Adair^a, Gordon Allan^a

^a Veterinary Sciences Division, Agri-Food and Biosciences Institute, Stormont, Belfast BT4 3SD, United Kingdom

^b Department of Veterinary Science, Queen's University of Belfast, Stormont, Belfast BT4 3SD, United Kingdom

^c Department of Virology, National Veterinary Institute, DTU, Lindholm, DK-4771 Kalvehave, Denmark

^d Centro de Investigación en Sanidad Animal (CISA-INIA), Valdeolmos, Madrid 28130, Spain

Article history: Received 8 October 2008 Received in revised form 30 April 2010 Accepted 6 May 2010 Available online 8 June 2010

Keywords: African swine fever virus Real-time PCR Minor groove binder probe

ABSTRACT

The design of a 5' conjugated minor groove binder (MGB) probe real-time PCR assay is described for the rapid, sensitive and specific detection of African swine fever virus (ASFV) DNA. The assay is designed against the 9GL region and is capable of detecting 20 copies of a DNA standard. It does not detect any of the other common swine DNA viruses tested in this study. The assay can detect ASFV DNA in a range of clinical samples. Sensitivity was equivalent to the Office International des Epizooties (OIE) recommended TaqMan assay. In addition the assay was found to have a detection limit 10-fold more sensitive than the conventional PCR recommended by the OIE. Linear range was ten logs from 2×10^1 to 2×10^{10} . The assay is rapid with an amplification time just over 2 h. The development of this assay provides a useful tool for the specific diagnosis of ASF in statutory or emergency testing programs or for the detection of ASFV DNA in research applications.

© 2010 Elsevier B.V. All rights reserved.

1. Introduction

African swine fever (ASF) is an OIE notifiable infectious disease of domestic and wild pigs and is caused by the African swine fever virus (ASFV), the sole member of the genus Asfivirus, family Asfarviridae. It is one of the most serious transboundary animal diseases and has the potential for rapid international spread with crippling socio-economic consequences. The disease is enzootic in sub-Saharan Africa, Madagascar and Sardinia where it can exist as a sylvatic cycle between warthogs and ticks of the genus Ornithodoros. The disease was first described in Kenya (Montgomery, 1921). The first European outbreak was in Lisbon in 1957 and since then outbreaks have occurred in several European countries, as well as in Cuba and Central and South America. In recent years it has spread through Georgia to several Caucasian countries. Slaughter of swine is the only means of eradication as there is no current treatment or vaccine available.

ASF is characterised by high fever of up to 42 °C, subcutaneous haemorrhage, inactivity, breathing difficulty, loss of appetite and often massive mortality. The mortality can vary as there is a range

of severity of disease, from chronic to peracute, depending on virus strain and host.

Deliberate or accidental feeding of pigs with human food waste is considered to be the most likely cause of introduction of the virus, particularly as ASFV can persist in dried meats for long periods (Mebus, 1988; McKercher et al., 1978). It is also known to survive for 18 months in proteinaceous substances at room temperature and for 11 days in non-proteinaceous substances such as faeces.

Diagnosis is complicated by the similarity of a range of other infections, particularly classical swine fever (CSF) (Kleiboeker, 2002) and porcine dermatitis and nephropathy syndrome (PDNS) (Chae, 2005). Rapid and reliable differential diagnosis is critical for the identification and control of the virus in the case of an outbreak and the OIE describes reliable and sensitive virus detection techniques for ASFV, including, real-time PCR and conventional PCR assays, in its Manual of Diagnostic Tests and Vaccines for Terrestrial Animals 2009, Chapter 2.8.1 (http://www.oie.int/).

Minor groove binder (MGB) probes (Afonina et al., 2002; Belousov et al., 2004) are short oligonucleotides with a minor groove binder molecule and a quencher conjugated to the 5' end, and a fluorophore to the 3' end. The minor groove binder molecule functions by intercalating into the minor groove of the DNA α -helix and stabilises the bond, thereby increasing the melting temperature (*Tm*) of the probe and allowing the use of shorter probes. The fluorophore and the quencher are held in close proximity to

^{*} Corresponding author. Tel.: +44 2890 525752; fax: +44 2890 525823. *E-mail address:* john.mckillen@afbini.gov.uk (J. McKillen).

^{0166-0934/\$ -} see front matter © 2010 Elsevier B.V. All rights reserved. doi:10.1016/j.jviromet.2010.05.005

each other due to the 3D conformation of the probe in solution. A fluorescent signal is triggered by hybridisation of the probe to a complementary target due to the resultant separation of fluorophore and quencher.

5' MGB probes have been used for the detection of hepatitis A virus, enteroviruses, influenza A and B, respiratory syncitial virus and monkeypox (Gersberg et al., 2006; Hymas and Hillyard, 2009; Li et al., 2006). This work describes the development of an MGB probe assay for the rapid and sensitive detection of ASFV. The 9GL region to which it has been designed is a promising target for ASFV detection due to the high levels of nucleotide conservation.

2. Materials and methods

2.1. Virus isolates and clinical samples

DNA from ASFV strain Spain 75 was used to develop the 9GL MGB assay. DNA extracts from 14 other isolates of ASFV were used to evaluate the assay. These isolates represent genogroups I, V, VIII and IX according to p72 subtyping (Gallardo et al., 2009) (Table 2).

Clinical samples originated from six pigs experimentally infected with the ASFV strain Malta/78 (pigs 1–3) or strain MwLil 20/1 (pigs 4–6). These isolates are ASFV p72 genogroups I and VIII respectively. Whole blood was taken from these animals at –1, 2, 5, 8, 9, 11 and 12 days post-infection. The pigs were euthanized between days 5 and 12 post-infection and samples of spleen and musculus quadriceps were taken from each. Meat juice was obtained by gentle freezing and thawing of the meat. Details of these samples are in Tables 4a and 4b.

DNA from a range of Aujeszky's disease virus (ADV), porcine circovirus type 2 (PCV2) and porcine parvovirus (PPV) isolates was selected for specificity testing. In addition, one bovine parvovirus type 1 (BPV1), porcine circovirus type 1 (PCV1), and one swine adenovirus isolate were selected. DNA from fetal calf lung and pig kidney (PK15A) cell lines was used for negative samples. Details of these samples are in Table 2.

2.2. Design of primers and probes

An alignment of 22 Genbank sequences (http://www.ncbi. nlm.nih.gov) from the 9GL region of the ASFV genome was created using the GeneDoc program (www.nrbsc.org/ gfx/genedoc/index.html). This area of the genome was chosen as it represented an area of high sequence homology between isolates. These sequences cover a geographical and chronological range of ASFV isolates. GenBank accession numbers and details of these sequences are shown in Table 1.

A number of MGB probes (Quantiprobes) and primer sets were designed using dedicated software (Qiagen, www.qiagen.com/). A suitable 17 mer MGB probe was chosen in an area of complete conservation on the alignment and a compatible primer set was chosen giving an amplicon size of 172 base pairs.

Probe sequence is 5'-CTGAAAGTCCTCCGAGT-3', forward primer sequence is 5'-GTTGTTATGGAACGCGAAG-3', and reverse primer sequence is 5'-CGCTCCTAGCTGGAAAGAAAA-3'.

A set of primers spanning the assay region to produce standard DNA for optimisation was designed using Oligo 5.0 (Molecular Biology Insights). Primer and probe sequences were subject to BLAST analysis to determine potential cross homologies with nucleotide sequences in the Genbank database (http://www.ncbi.nlm.nih.gov/). The MGB probe was purchased from Qiagen (Qiagen L.T.D., Crawley, U.K.) and all the primers used in the study were from Eurofins MWG Operon (Eurofins MWG Operon, London, U.K.).

2.3. Extraction of DNA

Nucleic acids from the 14 ASFV strains were extracted by Nucleospin kit (Clontech, Saint-Germain-en-Laye, France) according to the manufacturer's instructions.

DNA purification of ASFV from Spain 75, Malta and Malawi clinical samples in GuSCN was performed using neutral phenol to preclean the blood and meat juice samples, followed by Boom silica extraction of DNA as described for RNA extraction (Uttenthal et al., 2003). The DNA was eluted in water and adjusted to 0.3 M sodium acetate; 2 volumes of ethanol were added and the DNA was stored frozen for 24h before shipment to our laboratories. The samples were precipitated using a standard protocol. ADV, BPV1, PCV1, PCV2, porcine adenovirus and cell culture nucleic acids were extracted from 200 μ l amounts using Qiagen DNA Blood Mini Kit (Qiagen, Paisley, U.K.) according to the manufacturers' instructions.

Table 1

Details of GenBank ASFV sequences in the multiple sequence alignment used for MGB assay design.

Accession no.	Isolate	Year	Virulence	Origin	vp72 genotype
AF081173	Wildebeeslaagte/96/1 (M1)	1996	Virulent	South Africa	-
AF081163	Cr96/3	1996	Virulent	South Africa	-
AF081178	Ten62	1962	Virulent	Malawi	-
AF081176	Pretorisuskop/96/5 (Pr5)	1996	-	South Africa	-
AF081174	Mw Lil20/1ª	1983	Virulent	Malawi	VIII
AF081168	Killean III	1959	-	Kenya	-
AF081165	E75 ^a	1975	Virulent	Spain	Ι
NC_001659	BA71V ^a	1971	-	Spain	Ι
AF081160	Cameroon 82	1982	-	Cameroon	I
AF081170	Granja63	1963	Virulent	Spain	-
AF081172	Lis60 ^a	1960	Virulent	Portugal	I
AF081181	Kat67	1967	Virulent	Democratic Republic of the Congo	-
AF081169	Ken64	1964	Virulent	Kenya	-
AF081180	VF83	1983	Virulent	Zimbabwe	-
AF081171	Ken55	1955	Virulent	Kenya	-
AF081179	Ug61	1961	Virulent	Uganda	-
AF081177	Sp51	1951	Virulent	South Africa	-
AF081166	Haiti 81 ^a	1981	-	Haiti	I
AF081164	E70 ^a	1970	Virulent	Spain	I
AF081175	Pr4	1996	Virulent	South Africa	-
AF081162	Cr96/1	1996	Virulent	South Africa	-
AF081161	Ch83/1	1983	Virulent	Zimbabwe	-

^a Isolates available in this study.

2.4. Conventional gel-based PCR

Conventional gel-based PCR for the detection of ASFV nucleic acid was carried out using HotStarTaq Master Mix (Qiagen L.T.D., Crawley, U.K.) according to the manufacturer's instructions. The primers, used at a concentration of 0.5 µM, were those recommended by the OIE, Manual of Diagnostic Tests and Vaccines for Terrestrial Animals (ASFV OIE F 5'-ATGGATACCGAGGGAATAGC-3' and ASFV OIE R 5'-CTTACCGATGAAAATGATAC-3') and are designed against the p72 region of the ASFV genome. Thermal cycling was carried out on a DNA Engine Dyad thermal cycler (Bio-Rad Laboratories Ltd., Herts, U.K.). Cycling was initiated with a hot start denaturation of 15 min. Forty cycles followed consisting of 30 s each of denaturation at 94 °C, annealing at 50 °C, and extension at 72 °C. Reaction volumes were 25 µl. PCR products were visualized by electrophoresis on 1.5% (w/v) agarose gels (Mast Group L.T.D., Bootle, U.K.) and bands were photographed and annotated using BioRad GelDoc XR and Quantity One software (BioRad, Hemel Hempstead, U.K.).

2.5. Real-time MGB PCR

JumpStart Taq master mix (Sigma–Aldrich, Gillingham, U.K.) was used for the ASFV MGB assay. Probe concentration was $1 \times$ as recommended by the manufacturer (no details of actual concentrations are supplied). Primer concentrations of $0.5 \,\mu$ M each were used and Mg²⁺ was used at a final concentration of 3 mM. The assay was run on an MJ Research Opticon 2 real-time thermal cycler (Bio-Rad Laboratories Ltd., Herts, U.K.). The real-time PCR was initiated with a hot start denaturation of 2 min. This was followed by 40 cycles consisting of 30 s each of denaturation at 94 °C, annealing at 53 °C and extension at 72 °C. Reaction volumes were 25 μ l including 2 μ l of target. The MGB probe was FAM (3') labeled and quenched with an Eclipse Dark Quencher (5').

2.6. Real-time TaqMan PCR

Real-time TaqMan PCR was carried out using $2\times$ Gene Expression Mastermix (Applied Biosystems, Warrington, U.K.) according to the manufacturer's instructions. Primers at a final concentration of $2\,\mu$ M (5'-CTGCTCATGGTATCAATCTTATCGA, 5'-GATACCAACAAGATC(AG)GCCGT-3') and a TaqMan probe at a final concentration of $0.2\,\mu$ M (5'-[6-carboxy-fluorescein (FAM)]-CCACGGGAGGAATACCAACCAAGTG-3'-[6-carboxy-tetramethyl-rhodamine (TAMRA)]) were used. The reactions were made up to $23\,\mu$ l with sterile water and $2\,\mu$ l of target DNA was used as template. Cycling consisted of 95 °C for 10 min, and 40 cycles at 95 °C for 15 s and 60 °C for 1 min.

2.7. Production of dilution series of quantified standard DNA for assay optimisation

Bulk amounts of PCR product were used to produce a dilution series of quantified standard from strain Spain 75 DNA. A 433 base pair PCR amplicon was produced with a set of primers (F 5'-CCGGGAGACGTTGTTTT-3' and R 5'-TGTTGCATTGGGGACCTA-3') that spanned the assay region. PCR was carried out with an annealing temperature (*Ta*) of 52 °C and as described for the conventional ASFV PCR. The PCR product was electrophoresed on a 1.5% agarose gel (Mast Group L.T.D., Bootle, U.K.), the bands were excised, purified from the gel using QIAEX Gel Extraction Kit (Qiagen L.T.D., Crawley, U.K.) and quantified using a Jenway Genova UV/VIS spectrophotometer (Barloworld Scientific Ltd., Essex, U.K.). Based on the size of the amplicon and on the average molecular weight of a base pair the estimated number of copies of the amplicon/µl was calculated. From this a dilution series of known amplicon copy numbers from 10^{10} to 10^0 copies/µl was produced.

2.8. Optimisation of the MGB assay and determination of sensitivity and linear dynamic range

The first stage of optimisation involved performing a gradient PCR with *Ta* and Mg²⁺ concentrations between 50–60 °C and 3–7 mM respectively. On determining the optimal *Ta*/Mg²⁺ combination a primer matrix was carried out. This involved varying the forward and reverse primers in all combinations of 0.1, 0.25, 0.5, 0.75 and 1.0 μ M concentrations. The MGB probe is supplied by the manufacturer as a 20× solution. Concentrations of 0.25×, 0.5×, 1× and 2× were tested with the optimal *Ta*, Mg²⁺ and primer concentrations. Optimisation was carried out using both 2 × 10² and 2 × 10⁶ copies/ μ l of target standard DNA. Dynamic range and sensitivity of the optimised assays were determined by running triplicates of the quantified Spain 75 DNA standard material ranging from 2 × 10⁰ to 2 × 10¹⁰ copies/ μ l.

2.9. Specificity testing

To test specificity the ASFV MGB assay was run in parallel to the conventional PCR and TaqMan assays using nucleic acid from the ADV, ASFV, BPV1, PCV1, PCV2, porcine adenovirus and PPV virus isolates previously described (see Table 2). In addition pig kidney (PK15A) and fetal calf lung cell lines were used as controls.

2.10. Comparison of the MGB assay sensitivity to OIE recommended gel-based PCR and TaqMan PCR

In order to compare sensitivity of the MGB assay to gel-based and TaqMan PCR a DNA extract from ASFV isolate Haiti 1981 was diluted 10-fold in DEPC treated water (Ambion, Cambridgeshire, U.K.). Five replicates of each dilution were run with all three assays.

2.11. Comparison of MGB assay and OIE recommended gel-based PCR and TaqMan PCR for the detection of ASFV DNA in clinical samples

Whole bloods, spleen and meat juice samples were tested for the presence of ASFV DNA using the MGB assay and both the OIE recommended conventional PCR and TaqMan assays as described (Tables 4a and 4b).

3. Results

3.1. Optimisation of the MGB assay, determination of sensitivity and linear dynamic range

Optimisation of each assay involved determining the best Ta, Mg^{2+} concentration, primer concentrations and probe concentration. In a real-time PCR assay the cycle number at which the fluorescent trace for a particular reaction crosses a pre-determined threshold level is known as the Ct value. Conditions for the MGB assay were chosen that gave the lowest Ct values for a particular known starting concentration of template. In this way an assay of the greatest sensitivity was developed. The assay worked optimally with 3 mM Mg²⁺, a *Ta* of 53 °C, primer concentrations of 0.5 μ M each and a probe concentration of 1×. It was observed that the assay functioned well over a wide range of *Tas*. Cts for each temperature between 50 and 57 °C were within a one Ct range.

The assay was capable of detecting triplicates of this standard DNA across a ten orders of magnitude range, from 2×10^1

Table 2

Details of virus isolates tested in this study and results of specificity testing using ASFV MGB, TaqMan and conventional PCR assays. All ASFV accession numbers refer to vp72 gene.

ADVApy Total imp 76 imp 86 imp 87 imp 86 imp 86 imp 87 imp 87 imp 86 imp 87 imp 87 imp 87 imp 87 imp 88 imp 88 im		ISOLATE	Origin	Accession no.	p72 genotype	Ct	Ct	
ADVImp 76 Imp 7						ASFV MGB	ASFV TM	
Imp 55 Imp 56 Imp 5	ADV	Imp 76	England	-	-	No Ct	No Ct	-ve
Imp 156DemarkNo CtNo C		Imp 95	Holland	-	-	No Ct	No Ct	-ve
Imp 186 Imp 186 Imp 188 Imp 189 Imp 189 Im		Imp 156	Denmark	-	-	No Ct	No Ct	-ve
Inp 187BeigiumNo CtNo C		Imp 186	Belgium	-	-	No Ct	No Ct	-ve
Imp 188 Imp 209 Imp 209 BarhaBelgiumNo CiNo Ci-veBarha Burha Van DornHungaryNo CiNo CiNo Ci-veNa DornHollandNo CiNo Ci-ve-veASF Muli 2010HollandNo CiNo Ci-ve-veASF Muli 2011Moraben ireland Angola 1970F]174376I1131+veAngo Cape Verel 1970 CY070 Cape Verel 1970F]174380I2627+veAngola 1971 Chalswa 1983N261361Via2627+veVia DornChalswa 1983N216161Via2627+veCY070 Cape Verel 1970F]154413I2627+veKen07.Eld1Kenya 2005F]154413IX2627+veF70 E70Burkin Caso 2007F]154413IX2627+veF70 E70Burkin Caso 2007AY578692I242425+veF70 E70Barlinis I1975AY578692I2422+ve+veF70 E70Barlinis I1975AY578692I2422+ve+veF70 E70Barlinis I1975AY578692I2426-ve+veF71 E71Barlinis I1975AY578692I2426-ve+veF71 E71Barlinis I1971AY578693I		Imp 187	Belgium	-	-	No Ct	No Ct	-ve
Imp 209USÅNo CiNo Ci-ve-veBarhaHungaryNo CiNo CiNo Ci-ve-veNIA3Norther IrelandNo CiNo Ci-ve-veASFVMoraMorambique 1964F174376V3131-ve-veMuralAngola 1972F174378VII2624+ve-veMuli 2011Chalasva 1983AY261361VII2624+veMuli 2014Hoima, Uganda 2003F1154434IX3132+veKen07.Eid1Kenya 2007F1154434IX3132+veB670Porteverda Spain 1970F1154434IX3132+veB701Porteverda Spain 1970F1154434IX3132+veB670Porteverda Spain 1970PI154434IX3132+veB701Porteverda Spain 1970F1154345I2024+veB701Porteverda Spain 1970AY578693I2324+veB711Libon Portuga1060AY501530I2427+veB701Icition Spain 1975AY578693I2628+veB701Icition Spain 1975H373872I2628+veB701Icition Spain 1975AY578693I2628+veB701Icition Spain 1975AY578693IA		Imp 188	Belgium	-	-	No Ct	No Ct	-ve
Image: Section of the section of th		Imp 209	USA	-	-	No Ct	No Ct	-ve
IMA3 Na DoornNoft HollandNo Ct No CtNo Ct N		Bartha	Hungary	-	-	No Ct	No Ct	-ve
Na DoomHolland-No CtNo Ct-eASPMozé4Magola 1964F1/4376V3131**eAngola 1972F1/14376V1031**e**eMyUL20/1Chalswa 1983X20161VIII2627**eUG03HCape Verol 1997F1/5428VIII2627**eUg03HKeing 2006F1/54428IX2929**eKeing EllKeing 2007F1/54434IX3132**eBf07Burking Easo 2007F1/54434IX2125**eF07Ponteverde Spain 1970AV57869I2022**eBf07Burking Easo 2007F1/54434IX2022**eF07Burking Easo 2007F1/5438I2022**eF07Burking Easo 2007F1/5436F1/5438I2020**e		NIA3	Northern Ireland	-	-	No Ct	No Ct	-ve
ASFVMoz64 Ang72Mozambigue 1964 Ang72F174376 1174378V1111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111111		Van Doorn	Holland	-	-	No Ct	No Ct	-ve
Ang72Ang0a 1972[174378I2624+veMull 20/1Chalswa 1983X20161VII2627+ve(V97Cap Verde 1997[174380I2932+veUg03HHoima, Uganda 2003[154434IX3132-veKen06.B1Kenya 2006[154434]IX3132+veKen07.Eld1Kenya 2007-I2425+veB70Burkino Fao 2007-I2425+veF70Pontevedra Spain 1970AY578692I2324+veB71Badaic Spain 1971FI74348I2022+veE75Lérida Spain 1975AY578693I2425+veE76Badaic Spain 1975AY578693I2729+veE75Lérida Spain 1975AY578693I2729+veE76Badaic Spain 1975AY578693I2729+veF07Badaic Spain 1975AY578693I2729+veF08Sardina, Italy 1988FI74362I2628+veF07WeibridgeSardina, Italy 1988FI74362I2628+veF07DemarkAP655391-AF07SpainSardinaNo CtNo Ct-veF07SpainCanadaAP655391-No Ct <td>ASFV</td> <td>Moz64</td> <td>Mozambique 1964</td> <td>FJ174376</td> <td>V</td> <td>31</td> <td>31</td> <td>+ve</td>	ASFV	Moz64	Mozambique 1964	FJ174376	V	31	31	+ve
Image: here in the image: here in t		Ang72	Angola 1972	FJ174378	Ι	26	24	+ve
CV97Cope Verde 1997F]174380I2932+veU03HHoima, Uganda 2003F]154428IX3132+veKen06.B1Kenya 2007F]154441IX3132+veKen07.Eld1Kenya 2007F]154441IX2627+veF00Burkino Faso 2007-I242627+veF07Burkino Faso 2007AV578692I2324+veF07Burkino Faso 2007AV578693I2324+veF07Burkino Faso 2007AV578693I2427+veF07Badjao Spain 1971FJ74380I2022+veF07Badjao Spain 1975AV578693I2628+veF07Badjao Spain 1975AV578693I2628+veF07Badjao Spain 1975AV578693I2628+veF07Badjao Spain 1975AV578693I2628+veF07Badjao Spain 1975AV578693I2022+veF07Badjao Spain 1975AV578693I2627+veF07Badjao Spain 1975AV578693I2627+veF07Badjao Spain 1975AV578693I2626-veF07Badjao Spain 1975AV578693IA2022+veF07Badjao Spain 197FJ74350II <t< td=""><td></td><td>MwLil 20/1</td><td>Chalaswa 1983</td><td>AY261361</td><td>VIII</td><td>26</td><td>27</td><td>+ve</td></t<>		MwLil 20/1	Chalaswa 1983	AY261361	VIII	26	27	+ve
Image: big		CV97	Cape Verde 1997	FJ174380	Ι	29	32	+ve
keno6.81 keno7.2kl1Kenya 2006 kenya 2007 kenya 2007 kenya 2007Fj15434 Fj15444IX3132+veKeno7.2kl1Kenya 2007 kenya 2007F12627+veF70 F70 F70Dentevedra Spain 1970 R570AY578692I2324+veF70 F75 L60Leida Spain 1975 Lisbon Portugal 1960AF301530I2427+veF75 L60Lisbon Portugal 1960 Haiti 181AF301530I2427+veF70 L60Lisbon Portugal 1960AF301530I2628+veF70 Port-au-Prince 81Englan02022+veF701WeybridgeEnglan0No CtNo Ct-veF7021010 Stoon L1015 StoonCanadaAF055391 AF055391-No CtNo Ct-ve1013Spain DenmarkNo CtNo Ct-ve-ve1016SpainNo CtNo Ct-ve1017 L206BelgiumNo CtNo Ct-ve1018SpainNo CtNo Ct-ve1019SvedenNo CtNo Ct-ve1016BelgiumNo CtNo Ct-ve1016SvedenNo CtNo Ct-ve1017SvedenNo CtNo Ct-ve1018Belgiu		Ug03H	Hoima, Uganda 2003	FJ154428	IX	29	29	+ve
ken07.Eld1Kenya 2007Fj15441IX2627+ve +ve 1BF07Burkino Faso 200712425+ve +ve 1BF07Portevedra Spain 1970AY57869212427+ve 1Ba71VBadajoz Spain 1971FJ74378312022+ve 1Ba71VLérida Spain 1975AY57869212427+ve 1L60Lisbon Portugal 1960AF30153912729+ve 1S58Sardinia, Italy 1988FJ17436212628+ve 1PCV1WeytigeEngland2628+ve 1PCV299USAAF055391-No CtNo Ct-ve 1010 Stoon-ve 1017Pomark-1010 StoonCanadaAF055392-No CtNo Ct-ve 1017-ve 1016No Ct-ve 1017-ve 10191010 StoonEnglandNo CtNo Ct-ve 1016-ve 1016-ve 1016-ve 1017-ve 1016-ve 1016-ve 1016-ve 1017-ve 1016-ve 1016-ve 1016-ve 1017-ve 1016-ve 1016-ve 1016-ve 1017-ve 1016-ve 1016-ve 1016-ve 1016-ve 1016-ve 1016-ve 1016-ve 1016-ve 1016-ve 1016-ve 1016-ve 1016-ve 1016-ve 1016		Ken06.B1	Kenya 2006	FJ154434	IX	31	32	+ve
BF07Burkino Faso 2007-I2425+veF70Pontevedra Spain 1970P(174348)12022+veB71VBadajoz Spain 1975P(174348)12022+veF75Lérida Spain 1975AY57869212427+veE75Lérida Spain 1975AY57869212628+veF75Lérida Spain 1975P(17432)12022+veS88Sardinia, Italy 1988F(17432)12022+vePCV1WeybridgeEnglandNo CtNo Ct-vePCV2999USAAF055392-No CtNo Ct-ve1010 StoonCanadaAF055392-No CtNo Ct-ve1017DenmarkNo CtNo Ct-ve1016 StoonEnglandNo CtNo Ct-ve1016 StoonEnglandNo CtNo Ct-ve1016 StoonEnglandNo CtNo Ct-ve1026BelgiumNo CtNo Ct-ve1047SvedenNo CtNo Ct-ve1058FranceAF055394-No CtNo Ct-ve1047SvedenNo CtNo Ct-ve1058FranceNo CtNo Ct-ve1059Fran		Ken07.Eld1	Kenya 2007	FJ154441	IX	26	27	+ve
E70Pontevedra Spain 1970AYS7892I2324+veBa71VBadajoz Spain 1975F]174342I2092+veBA75Lerida Spain 1975AYS78693I242792+veL60Lisbon Portugal 1960AF301539I2628+vePOrta-u-Prince 81Hit 1981F]174362I202020-vePCV1WeybridgeEnglandNo CtNo Ct-vePCV299USAAF055391-No CtNo Ct-ve1010 StoonCanadaAF055392-No CtNo Ct-ve1019DenmarkNo CtNo Ct-ve1019DenmarkNo CtNo Ct-ve1019DenmarkNo CtNo Ct-ve1019SpainNo CtNo Ct-ve1019SpainNo CtNo Ct-ve1026EnglandNo CtNo Ct-ve1036EnglandNo CtNo Ct-ve1247SwedenNo CtNo Ct-ve1247SwedenNo CtNo Ct-ve1247SwedenNo CtNo Ct-ve1247SwedenNo CtNo Ct-ve1247Sweden<		BF07	Burkino Faso 2007	_	Ι	24	25	+ve
Ba71VBadajoz Spain 1971FIJ 74348I2022+veE75Lérida Spain 1975AY5809I2427+veE60Lison Portugal 1960AF301539I2729+veSx80Sardinia, Italy 1988FIJ 74362I2628+vePCV1WeybridgeEnglandRo CtNo Ct-vePCV2999USAAF055391-No CtNo Ct-ve1010 StoonCanadaAF055391-No CtNo Ct-ve1017DenmarkNo CtNo Ct-ve1019SpainNo CtNo Ct-ve1019SpainNo CtNo Ct-ve1036EnglandNo CtNo Ct-ve1046EnglandNo CtNo Ct-ve1056EnglandNo CtNo Ct-ve1047SwedenNo CtNo Ct-ve1048HaugaryNo CtNo Ct-ve1049No CtNo Ct-ve1049No CtNo Ct-ve1050HagaryNo CtNo Ct-ve1051HagaryNo Ct-ve1051EnglandNo Ct-ve-ve1		E70	Pontevedra Spain 1970	AY578692	Ι	23	24	+ve
F75Lérida Spain 1975M578633I2427+veL60Lisbon Portugal 1960AF301539I2728*veS588S588T174362I2628*vePCV1WeybridgeEnglandNo CtNo Ct-ePCV2999USAAF055391-No CtNo Ct-ve1010CanadaAF055392-No CtNo Ct-ve1017DenmarkNo CtNo Ct-ve1018SpainNo CtNo Ct-ve1019SpainNo Ct-ve1019SpainNo Ct-ve1016SpainNo Ct-ve1018SpainNo Ct-ve1019SpainNo Ct-ve1018SpainNo Ct-ve1019SpainNo Ct-ve1026BelgiumNo Ct-ve1030TalgandNo Ct-ve1247SwedenNo CtNo Ct-ve1247SwedenNo CtNo Ct-ve1247SassFranceNo Ct-ve1247SassFranceNo Ct-ve1247SassFranceNo		Ba71V	Badajoz Spain 1971	FJ174348	Ι	20	22	+ve
Loo Ss8 Ss8 Pr1-au-Prince 81Lisbon Portugal 1960 Sardinia, Italy 1988AF301539 F1743C2 F1743C2I27 26 26 27 2029 20+ve 20PCV1WeybridgeEngland2022+ve +vePCV299.9USAAF055391-No CtNo Ct-ve -ve1010 StoonCanadaAF055392-No CtNo Ct-ve -ve1017DenmarkNo CtNo Ct-ve -ve1019SpainNo CtNo Ct-ve -ve1019SpainNo CtNo Ct-ve -ve1019SpainNo CtNo Ct-ve -ve1019SpainNo CtNo Ct-ve -ve1019SpainNo CtNo Ct-ve -ve1019SpainNo CtNo Ct-ve -ve1026EnglandNo CtNo Ct-ve -ve1036EnglandNo CtNo Ct-ve -ve1042FranceNo CtNo Ct-ve -ve1057FranceNo CtNo Ct-ve -ve105SpafFranceNo CtNo Ct-ve -ve105SpafNo CtNo Ct-ve -ve105Spaf<		E75	Lérida Spain 1975	AY578693	Ι	24	27	+ve
S888 Port-au-Prince 81Sardinia, Ital y 1988 Hait i 1981FJ174322 FJ174375I26 2028 22+ve +vePCV1WeybridgeEnglandNo CtNo Ct-vePCV2999USAAF055391-No CtNo Ct-ve1010 StoonCanadaAF055392-No CtNo Ct-ve1017DenmarkNo CtNo Ct-ve1019SpainNo CtNo Ct-ve1206BelgiumEF990644-No CtNo Ct-ve1206BelgiumEF990644-No CtNo Ct-ve1206BelgiumNo CtNo Ct-ve1247SwedenNo CtNo Ct-ve1247SwedenNo CtNo Ct-ve9367ItalyNo CtNo Ct-ve9367ItalyNo CtNo Ct-ve9367ItalyNo CtNo Ct-ve9367ItalyNo CtNo Ct-ve1005NDLSandaNo CtNo Ct-ve1005ShapeFranceNo CtNo Ct-ve1007SeeSeeNo CtNo Ct-ve1008SeeNo CtNo CtNo Ct-ve-ve		L60	Lisbon Portugal 1960	AF301539	Ι	27	29	+ve
Port-au-Prince 81Haiti 1981FJ174375I2022+vePCV1WeybridgeEnglandNo CtNo Ct-vePCV2999USAAF055391-No CtNo Ct-ve1010 StoonCanadaAF055392-No CtNo Ct-ve1017DenmarkNo CtNo Ct-ve1019SpainNo CtNo Ct-ve1206BelgiumEF990644-No CtNo Ct-ve1306EnglandNo CtNo Ct-ve1247SwedenNo CtNo Ct-ve5549Northern IrelandNo CtNo Ct-ve9367ItalyNo CtNo Ct-ve9367HugaryNo CtNo Ct-ve1008FranceAF055394No CtNo Ct-ve9367HugaryNo CtNo Ct-ve1008FranceNo CtNo Ct-ve1009SpainNo CtNo Ct-ve1009FranceNo CtNo Ct-ve1009SpainNo CtNo Ct-ve1009FranceNo CtNo Ct-ve1009FranceNo Ct		Ss88	Sardinia, Italy 1988	FJ174362	Ι	26	28	+ve
PCV1WeybridgeEnglandNo CtNo Ct-vePCV2999USAAF055391-No CtNo CtNo Ct-ve1010 StoonCanadaAF055392-No CtNo CtNo Ct-ve1017DenmarkNo CtNo CtNo Ct-ve1019SpainNo CtNo Ct-ve1009BelgiumEP90644-No CtNo Ct-ve1306EnglandNo CtNo Ct-ve1247SwedenNo CtNo Ct-ve1247SwedenNo CtNo Ct-ve9367ItalyNo CtNo Ct-ve9367ItalyNo CtNo Ct-ve12755FranceAF055394-No CtNo Ct-ve1005 NADLSpandaNo Ct-ve1005 NADLCanadaNo Ct-ve1005 NADLSpandave1005 NADLSpandaNo CtNo Ct-veSwine ademField caseNorter IrelandNo CtNo Ct-veSwine ademField caseNo CtSWINE ademNo CtNo Ct-veS		Port-au-Prince 81	Haiti 1981	FJ174375	Ι	20	22	+ve
PCV2999USAAF055391-No CtNo Ct-ve1010 StoonCanadaAF055392-No CtNo Ct-ve1017DemarkNo CtNo Ct-ve1019SpainNo CtNo Ct-ve1206BelgiumEF990644-No CtNo Ct-ve1306EnglandNo CtNo Ct-ve1247SwedenNo CtNo Ct-ve9367ItalyNo CtNo Ct-ve9367ItalyNo CtNo Ct-ve42285FranceAF055394-No CtNo Ct-ve1005 NADLEnglandNo CtNo Ct-ve1005 NADLEnglandNo CtNo Ct-ve1005 NADLEnglandNo CtNo Ct-ve1008FranceNo CtNo Ct-veSwine ademoField caseNorten IrelandNo CtNo Ct-veBPV1HadenNo CtNo Ct-ve-Swine ademoField caseNorten IrelandNo CtNo Ct-veFV14HadenNo CtNo Ct-ve-FV15HadenNo CtNo Ct-ve- <td>PCV1</td> <td>Weybridge</td> <td>England</td> <td>-</td> <td>-</td> <td>No Ct</td> <td>No Ct</td> <td>-ve</td>	PCV1	Weybridge	England	-	-	No Ct	No Ct	-ve
1010 Stoon Canada AF055392 - No Ct No Ct ve 1017 Denmark - - No Ct No Ct ve 1019 Spain - - No Ct No Ct ve 1019 Spain EF990644 - No Ct No Ct ve 1206 Belgium EF990644 - No Ct No Ct ve 1306 England - - No Ct No Ct ve 1247 Sweden - - No Ct No Ct ve 549 Northern Ireland - - No Ct No Ct -ve 9367 Italy - - No Ct No Ct -ve H2755 Hungary - - No Ct No Ct -ve PV 59 England - - No Ct No Ct -ve 1005 NADL Canada NC001718 -	PCV2	999	USA	AF055391	-	No Ct	No Ct	-ve
1017 Denmark - - No Ct No Ct -ve 1019 Spain - - No Ct No Ct -ve 1206 Belgium EF990644 - No Ct No Ct -ve 1306 England - - No Ct No Ct -ve 1247 Sweden - - No Ct No Ct -ve 5549 Northern Ireland - - No Ct No Ct -ve 9367 Italy - - No Ct No Ct -ve 48285 France AF055394 - No Ct No Ct -ve PPV 59e England - - No Ct No Ct -ve 1005 NADL Canada NC.001718 - No Ct No Ct -ve Swine adeno Field case Northern Ireland - - No Ct No Ct -ve PV1 Haden - <td></td> <td>1010 Stoon</td> <td>Canada</td> <td>AF055392</td> <td>-</td> <td>No Ct</td> <td>No Ct</td> <td>-ve</td>		1010 Stoon	Canada	AF055392	-	No Ct	No Ct	-ve
Independence Spain - - No Ct No Ct -ve 1206 Belgium EF990644 - No Ct No Ct -ve 1306 England - - - No Ct No Ct -ve 1306 England - - - No Ct No Ct -ve 1247 Sweden - - No Ct No Ct -ve 5549 Northern Ireland - - No Ct No Ct -ve 9367 Italy - - No Ct No Ct -ve 48285 France AF055394 - No Ct No Ct -ve PPV S9e England - - No Ct No Ct -ve 1005 NADL Canada NC.001718 - No Ct No Ct -ve Swine adeno Field case Northern Ireland - - No Ct No Ct -ve BP		1017	Denmark	-	-	No Ct	No Ct	-ve
1206 Belgium EF990644 - No Ct No Ct -ve 1306 England - - No Ct No Ct -ve 1247 Sweden - - No Ct No Ct -ve 5549 Northern Ireland - - No Ct No Ct -ve 9367 Italy - - No Ct No Ct -ve 48285 France AF055394 - No Ct No Ct -ve PPV 59e England - - No Ct No Ct -ve 1005 NADL Canada NC.001718 - No Ct No Ct -ve Swine adeno Field case Northern Ireland - - No Ct No Ct -ve BPV1 Haden - - No Ct No Ct -ve FSwine adeno Field case Northern Ireland - - No Ct No Ct -ve BPV		1019	Spain	-	-	No Ct	No Ct	-ve
1306 England - - No Ct No Ct -ve 1247 Sweden - - No Ct No Ct -ve 5549 Northern Ireland - - No Ct No Ct -ve 9367 Italy - - No Ct No Ct -ve 48285 France AF055394 - No Ct No Ct -ve 1005 NADL England - - No Ct No Ct -ve 1005 NADL Canada NC.001718 - No Ct No Ct -ve Swine adeno Field case Northern Ireland - - No Ct No Ct -ve PV1 Haden - - No Ct No Ct -ve Swine adeno Field case Northern Ireland - - No Ct No Ct -ve PK15A - - - - No Ct No Ct -ve FCL		1206	Belgium	EF990644	-	No Ct	No Ct	-ve
1247 Sweden - - No Ct No Ct -ve 5549 Northern Ireland - - No Ct No Ct -ve 9367 Italy - - No Ct No Ct -ve 48285 France AF055394 - No Ct No Ct -ve H2755 Hungary - - No Ct No Ct -ve PPV 59e England - - No Ct No Ct -ve 1005 NADL Canada NC.001718 - No Ct No Ct -ve Swine adeno Field case Northern Ireland - - No Ct No Ct -ve BPV1 Haden - - - No Ct No Ct -ve PK15A - - - - No Ct No Ct -ve FCL - - - - No Ct No Ct -ve		1306	England	-	-	No Ct	No Ct	-ve
5549Northern IrelandNo CtNo Ct-ve9367ItalyNo CtNo Ct-ve48285FranceAF055394-No CtNo Ct-veH2755HungaryNo CtNo Ct-vePPV59eEnglandNo CtNo Ct-ve1005 NADLCanadaNC_001718-No CtNo Ct-ve1008FranceNo CtNo Ct-veSwine adenoField caseNorthern IrelandNo CtNo Ct-veBPV1HadenNo CtNo Ct-vePK15ANo CtNo Ct-veFCLNo CtNo Ct-veFCL		1247	Sweden	-	-	No Ct	No Ct	-ve
9367ItalyNo CtNo Ct-ve48285FranceAF055394-No CtNo Ct-veH2755HungaryNo CtNo Ct-vePPV59eEnglandNo CtNo Ct-ve1005 NADLCanadaNC.001718-No CtNo Ct-ve1008FranceNo CtNo Ct-veSwine adenoField caseNorthern IrelandNo CtNo Ct-veBPV1HadenNo CtNo Ct-vePK15ANo CtNo Ct-veFCLNo CtNo Ct-ve		5549	Northern Ireland	-	-	No Ct	No Ct	-ve
48285 H2755FranceAF055394-No CtNo Ct-vePPV59e 1005 NADL 1008EnglandNo CtNo Ct-veFranceCanadaNC.001718No CtNo Ct-veSwine adenoField caseNorthern IrelandNo CtNo Ct-veBPV1HadenNo CtNo Ct-vePK15ANo CtNo Ct-veFCLNo CtNo Ct-ve		9367	Italy	-	-	No Ct	No Ct	-ve
PPVFugaryNo CtNo Ct-ve1005 NADL 1005 NADL 1008England CanadaNo CtNo Ct-veFranceNo CtNo Ct-veSwine adeno BPV1Field caseNorthern IrelandNo CtNo Ct-veSWINE adeno PV1No CtNo Ct-veBPV1 PK15A FCLNo CtNo Ct-veFCLNo CtNo Ct-veFCLNo CtNo Ct-ve		48285	France	AF055394	-	No Ct	No Ct	-ve
PPV59eEnglandNo CtNo Ct-ve1005 NADL 1008CanadaNC.001718-No CtNo Ct-ve1008FranceNo CtNo Ct-veSwine adenoField caseNorthern IrelandNo CtNo Ct-veBPV1HadenNo CtNo Ct-vePK15ANo CtNo Ct-veFCLNo CtNo Ct-ve		H2755	Hungary	-	-	No Ct	No Ct	-ve
1005 NADL 1008CanadaNC_001718-No CtNo Ct-ve1008FranceNo CtNo Ct-veSwine adenoField caseNorthern IrelandNo CtNo Ct-veBPV1HadenNo CtNo Ct-vePK15ANo CtNo Ct-veFCLNo CtNo Ct-ve	PPV	59e	England	-	-	No Ct	No Ct	-ve
1008FranceNo CtNo Ct-veSwine adenoField caseNorthern IrelandNo CtNo Ct-veBPV1HadenNo CtNo Ct-vePK15ANo CtNo Ct-veFCLNo CtNo Ct-ve		1005 NADL	Canada	NC_001718	-	No Ct	No Ct	-ve
Swine adenoField caseNorthern IrelandNo CtNo Ct-veBPV1HadenNo CtNo Ct-vePK15ANo CtNo Ct-veFCLNo CtNo Ct-ve		1008	France	-	-	No Ct	No Ct	-ve
BPV1 Haden - - No Ct No Ct -ve PK15A - - - No Ct No Ct -ve FCL - - - No Ct No Ct -ve	Swine adeno	Field case	Northern Ireland	-	-	No Ct	No Ct	-ve
PK15A - - No Ct No Ct -ve FCL - - - No Ct No Ct -ve	BPV1	Haden	_	-	-	No Ct	No Ct	-ve
FCL – – – No Ct No Ct -ve	PK15A	-	_	-	-	No Ct	No Ct	-ve
	FCL	-	-	-	-	No Ct	No Ct	-ve

to 2×10^{10} copies/µl, and showed good linearity in this range (Fig. 1).

3.2. Specificity testing

The MGB assay was tested using ADV, BPV1, ASFV, PCV2, PCV1, PPV and porcine adenovirus as well as PK15A and FCL cell line nucleic acids as target and only produced a fluorescent signal with its own specific target (Table 2).

3.3. Comparison of the MGB assay sensitivity to OIE recommended gel-based PCR and TaqMan PCR

Table 3 shows a side-by-side comparison of the three assays on a serial dilution of Haiti ASFV DNA. The two real-time assays showed equivalent sensitivity being able to detect 10^{-4} dilutions in 100% of the five replicates. The MGB assay detected 40% of the five replicates at 10^{-5} , whereas the TaqMan assay detected 20% of the replicates at this dilution. The conventional assay was one log less sensitive, being capable of detecting 100% of the replicates at a 10^{-3} dilution and 20% of the replicates at a 10^{-4} dilution. 3.4. Comparison of MGB assay and OIE recommended gel-based PCR and TaqMan PCR for the detection of ASFV DNA in clinical samples

Data was generated by the ASFV MGB assay run using DNA extracted from whole blood, meat juice and spleen from experimentally infected animals as described. This was compared to data generated using the OIE TaqMan and conventional PCR assays. Tables 4a and 4b show the results. The results from the three assays were extremely consistent. All whole blood samples that were

Table 3

Comparison of the sensitivity of the ASFV MGB, TaqMan and conventional PCR assays using titrated Haiti strain DNA as target.

Dilution	ASFV MGB	ASFV TaqMan	Electrophoresis
Haiti 10 ⁻²	+++++	+++++	+++++
Haiti 10 ⁻³	+++++	+++++	+++++
Haiti 10 ⁻⁴	+++++	+++++	+
Haiti 10 ⁻⁵	++	+	
Haiti 10 ⁻⁶			
Haiti 10 ⁻⁷			



Fig. 1. (a) Fluorescence vs. cycle no. for MGB assay run with triplicates of 2×10^{10} to 2×10^1 dilutions of standard DNA. (b) Slope and r^2 values from log quantity vs. Ct standard curve for the ASFV MGB assay with dilutions of standard DNA from 2×10^{10} to 2×10^1 , y = -0.28x + 11.58, and $r^2 = 0.99$.

detected by the TaqMan and conventional assays were detected by the MGB assay. Results for the spleen and meat juice samples were also largely consistent across all three assay formats. One spleen sample was negative for the conventional assay (pig no. 3, 11 days post infection) but positive for the two real-time assays.

4. Discussion

African swine fever is considered by the United Nations Food and Agriculture Organisation (FAO) to be a transboundary animal disease of high socio-economic importance based on its mortality, rapid international spread and the according threat to food and trade security. As ASFV is highly contagious and spreads rapidly, with a large amount of shedding (even before the onset of clinical signs) and up to 100% morbidity, it is highly important in the case of an outbreak for laboratories to confirm diagnosis as soon as possible. As such the FAO includes ASF in the emergency prevention of transboundary animal and plant pests and diseases (EMPRES). The OIE classes ASF as a notifiable disease. The OIE prescribed test for international transfer is by enzyme-linked immunosorbant assay (ELISA). Diagnosis can also be by virus isolation (VI) in pig monocytes, fluorescent antibody testing or immunoblotting (OIE Manual of Diagnostic Tests and Vaccines for Terrestrial Animals 2009, Chapter 2.8.1). Several conventional PCR assays have been developed for ASFV detection (Agüero et al., 2003; Steiger et al., 1992) including detection in Ornithodoros ticks (Basto et al., 2006). The OIE man-

Table 4a

Comparison of PCR results on DNA extracted from whole blood samples at various days post-infection (PI) from pigs infected with ASFV Malta/78 (pigs 1–3) or MwLil 20/1 (pigs 4–6) strains.

Pig no.	Days PI	MGB assay Ct	TaqMan assay Ct	OIE PCR
1	-1 2 5 8 12	No Ct No Ct No Ct No Ct 23	No Ct No Ct No Ct No Ct 23	Neg Neg Neg +
2	-1 2 5 8 9	No Ct No Ct 23 29 31	No Ct No Ct 27 33 34	Neg Neg + +
3	-1 2 5 8	No Ct No Ct 24 28	No Ct No Ct 25 29	Neg Neg + +
Pig no.	Days PI	MGB assay Ct	TM assay Ct	OIE PCR
Pig no. 3	Days PI 11	MGB assay Ct 37	TM assay Ct 37	OIE PCR +
Pig no. 3 4	Days PI 11 -1 2 5 8 9	MGB assay Ct 37 No Ct No Ct No Ct 29 25	TM assay Ct 37 No Ct No Ct No Ct 30 26	OIE PCR + Neg Neg Neg + +
Pig no. 3 4 5	Days PI 11 -1 2 5 8 9 -1 2 5	MGB assay Ct 37 No Ct No Ct 29 25 No Ct 32 22	TM assay Ct 37 No Ct No Ct 30 26 No Ct 32 22	OIE PCR + Neg Neg + + Neg + +

ual describes conventional and TaqMan ASFV PCR assays (King et al., 2003). Other TaqMan (Zsak et al., 2005) and Molecular Beacon (McKillen et al., 2007) methods have been described. The Invader assay, a novel isothermal amplification method using a cleavase enzyme has been applied to ASFV (Hjertner et al., 2005).

Real-time PCR and isothermal assays are now mainstream techniques for the detection of virus nucleic acids. Their speed, sensitivity and specificity have led to their widespread use for the detection of economically important viruses, including those of swine (Belák, 2005). The closed-tube format of real-time PCR assays such as this MGB assay means that no agarose gel electrophoresis is necessary, run times are shorter, potential for cross contamination is reduced and potential exposure to carcinogenic ethidium bromide is eliminated.

Accurate differential diagnosis of ASF is critical as it may be confused clinically with several other diseases particularly classical swine fever (CSFV) but also bacterial septicemias, fungal poisonings and viral diseases such as porcine reproductive and respiratory syndrome, porcine dermatitis and nephropathy syndrome and Aujeszky's disease (Kleiboeker, 2002). In the event of an outbreak even short delays in determining the correct diagnosis and

Table 4b

Comparison of PCR results on DNA extracted from spleen and meat juice samples from the pigs infected with ASFV Malta/78 (pigs 1–3) or MwLil 20/1 (pigs 4–6) strains.

Pig no.	Days PI	Spleen			Meat juice		
		MGB assay Ct	TaqMan assay Ct	OIE PCR	MGB assay Ct	TaqMan assay Ct	OIE PCR
1	12	24	26	+	32	36	+
2	9	24	25	+	34	33	+
3	11	36	35	-ve	33	34	+
4	9	23	23	+	33	34	+
5	5	19	19	+	27	27	+
6	8	20	21	+	29	29	+

initiating the appropriate response can have devastating effects on the abilities of agricultural authorities to curb the spread of infection. The MGB assay described here has a run time of just over 2 h. The speed of the assay is comparable to other real-time chemistries such as molecular beacons or TaqMan probes.

The 5' conjugated MGB probes (Afonina et al., 2002; Belousov et al., 2004) used in this study are short oligonucleotides with the MGB moiety attached to the 5' terminus. These are not to be confused with 3' conjugated TagMan-MGB probes, which are unconjugated at the 5' end to allow exonuclease digestion by the polymerase enzyme. The MGB molecule intercalates into the minor groove of the DNA double helix increasing the *Tm* of the complex by up to 30 °C and thereby allowing the design of shorter probes. This could be particularly useful in the design of assays for viruses that have considerable sequence divergence between strains such as many RNA viruses. The marked reduction in Tm of the probe/target hybrid when binding is to a mismatched target also confers a high level of specificity on the assay. These MGB probes can also utilise a number of modified bases that confer on them further diagnostic advantages. Both Super ATM and Super TTM, when incorporated into AT rich regions cause increased stabilisation. Super GTM bases allow the design of probes in high GC rich sequences. Normally design in such regions would be impossible due to the formation of tetrad aggregates of guanosine causing G-G association and an according reduction in hybridisation efficiency. The conserved nature of the ASFV genome in the 9GL region meant that incorporation of these modified bases was unnecessary for this assay. However, future work will include the evaluation of these probes for RNA viruses

The ASF MGB assay combines all the necessary characteristics for a modern molecular diagnostic test. The assay is highly sensitive, being capable of detection of as little as 20 copies of standard DNA per reaction. It can detect virus in whole blood, spleen and meat juice. All clinical samples detected by the conventional and TagMan OIE recommended assays were detected by this test. One of the whole blood sample that was detected by both the TagMan and the MGB assays was not detected by conventional PCR. This sample was weak (Ct values of 35-37) when run with the MGB and TagMan assays; so reduced sensitivity of the conventional PCR could explain this. Sensitivity testing showed that the two realtime assays were equally sensitive and that these could consistently detect ASFV Haiti to a dilution of 10⁻⁴. The conventional OIE assay was one 10-fold dilution less sensitive. Specificity testing on a range of DNA viruses showed that the MGB assay would only detect specific ASFV DNA and did not detect DNAs from ADV, BPV1, PCV2, PPV and swine adenovirus or from bovine and porcine cell lines. The ideal assay Ta was determined to be 53 °C. Optimisation showed little significant difference in functionality between Tas of 50 and 56 °C. As such for the detection of highly variable viral genomes a balance could be found through careful selection of Ta that allows tolerance of some mismatches whilst still retaining high specificity.

In many diagnostic situations high throughput pathogen detection is important. Real-time PCR methods allow such high throughput application. Real-time assays can be carried out in 96 well plate formats and this allows them to be adapted easily to automated platforms. Automation of both reaction set-up and nucleic acid extraction now allows for highly reproducible high throughput pathogen detection. MGB assays could readily be adapted for such automated processes. Cost of the assay probe is equivalent to other real-time reporter chemistries.

In summary we have developed an assay suitable for the rapid, reproducible and sensitive differential detection of ASFV DNA using MGB technology. This assay is suitable for applications in research into ASFV and could be used in emergency outbreak situations if validated further. The MGB technology was easily applicable to the molecular detection of ASFV and has unique qualities that could lend themselves to improved detection of more difficult viral targets.

Acknowledgements

This research formed part of an EU funded project entitled "New and Emerging Technologies: Improved Laboratory and On-Site Detection of OIE List A Viruses in Animals and Animal Products", *Website*: http://www.labonsite.com/. This collaborative project involved nine European laboratories and was supported by the 6th Framework Programme of the European Commission SSPE-CT-2004-513645.

The authors would like to thank Heidi Elbrink and Nethe Wrang Larsen, Department of Virology, National Veterinary Institute, Lindholm, Kalvehave, Denmark for excellent technical support.

References

- Afonina, I.A., Reed, M.W., Lusby, E., Shishkina, I.G., Belousov, Y.S., 2002. Minor groove binder-conjugated DNA probes for quantitative DNA detection by hybridizationtriggered fluorescence. Biotechniques 32 (4), 940–944. 946–949.
- Agüero, M., Fernández, J., Romero, L., Sanchez, C., Arias, M., Sánchez-Vizcaíno, J.M., 2003. Highly sensitive PCR assay for routine diagnosis of African swine fever virus in clinical samples. J. Clin. Microbiol. 41 (9), 4431–4434.
- Basto, A.P., Portugal, R.S., Nix, R.J., Cartaxeiro, C., Boinas, F., Dixon, L.K., Leitão, A., Martins, C., 2006. Development of a nested PCR and its internal control for the detection of African swine fever virus (ASFV) in *Ornithodoros erraticus*. Arch. Virol. 151, 819–826.
- Belák, S., 2005. The molecular diagnosis of porcine viral diseases: a review. Acta Vet. Hung. 53 (1), 113–124.
- Belousov, Y.S., Welch, R.A., Sanders, S., Mills, A., Kulchenko, A., Dempcy, R., Afonina, I.A., Walburger, D.K., Glaser, C.L., Yadavalli, S., Vermeulen, N.M., Mahoney, W., 2004. Single nucleotide polymorphism genotyping by two colour melting curve analysis using the MGB Eclipse Probe System in challenging sequence environment. Hum. Genomics 1 (March (3)), 209–217.
- Chae, C., 2005. A review of porcine circovirus 2-associated syndromes and diseases. Vet. J. 169 (May (3)), 326–336, Review.
- Gallardo, C., Mwaengo, D.M., Macharia, J.M., Arias, M., Taracha, E.A., Soler, A., Okoth, E., Martín, E., Kasiti, J., Bishop, R.P., 2009. Enhanced discrimination of African swine fever virus isolates through nucleotide sequencing of the p54, p72, and pB602L (CVR) genes. Virus Genes 38 (February (1)), 85–95.
- Gersberg, R.M., Rose, M.A., Robles-Sikisaka, R., Dhar, A.K., 2006. Quantitative detection of hepatitis a virus and enteroviruses near the United States–Mexico border and correlation with levels of fecal indicator bacteria. Appl. Environ. Microbiol. 72 (December (12)), 7438–7444.
- Hjertner, B., Meehan, B., McKillen, J., McNeilly, F., Belák, S., 2005. Adaptation of an invader assay for the detection of African swine fever virus DNA. J. Virol. Methods 124, 1–10.
- Hymas, W.C., Hillyard, D.R., 2009. Evaluation of nanogen MGB alert detection reagents in a multiplex real-time PCR for influenza virus types A and B and respiratory syncytial virus. J. Virol. Methods 156 (March (1–2)), 124–128.
- Kleiboeker, S.B., 2002. Swine fever: classical swine fever and African swine fever. Vet. Clin. North Am. Food Anim. Pract. 18, 431–451.
- King, D.P., Reid, S.M., Hutchings, G.H., Grierson, S.S., Wilkinson, P.J., Dixon, L.K., Bastos, A.D., Drew, T.W., 2003. Development of a TaqMan PCR assay with internal amplification control for the detection of African swine fever virus. J. Virol. Methods 107, 53–61.
- Li, Y., Olson, V.A., Laue, T., Laker, M.T., Damon, I.K., 2006. Detection of monkeypox virus with real-time PCR assays. J. Clin. Virol. 36 (July (3)), 194–203.
- McKercher, P.D., Hess, W.R., Hamdy, F., 1978. Residual viruses in pork products. Appl. Environ. Microbiol. 35 (1), 142–145.
- McKillen, J., Hjertner, B., Millar, A., McNeilly, F., Belák, S., Adair, B., Allan, G., 2007. Molecular beacon real-time PCR detection of swine viruses. J. Virol. Methods 40 (1–2), 155–165.
- Mebus, C.A., 1988. African swine fever. Adv. Virus Res. 35, 251-269, Review.
- Montgomery, R.E., 1921. On a farm of swine fever occurring in British East Africa (Kenya colony). J. Comp. Pathol. Ther. 34, 159–191, 243–264.
- Office International des Epizooties, 2009. Manual of diagnostic tests and vaccines for terrestrial animal. http://www.oie.int/.
- Steiger, Y., Ackermann, M., Mettraux, C., Kihm, U., 1992. Rapid and biologically safe diagnosis of African swine fever infection by using polymerase chain reaction. J. Clin. Microbiol. 30, 1–8.
- Uttenthal, Å., Storgaard, T., Oleksiewicz, M.B., de Stricker, K., 2003. Experimental infection with the Paderborn isolate of classical swine fever virus in 10-weekold pigs: determination of viral replication kinetics by quantitative RT-PCR, virus isolation and antigen ELISA. Vet. Microbiol. 92 (3), 197–212.
- Zsak, L., Borca, M.V., Risatti, G.R., Zsak, A., French, R.A., Lu, Z., Kutish, G.F., Neilan, J.G., Callahan, J.D., Nelson, W.M., Rock, D.L., 2005. Preclinical diagnosis of African swine fever in contact-exposed swine by a real-time PCR assay. J. Clin. Microbiol. 43, 112–119.